

JASCO

FVS-6000

Vibrational Circular Dichroism Spectrometer An essential tool for absolute configuration of molecules

JASCO has developed the FVS-6000 high performance VCD (Vibrational Circular Dichroism) spectrometer with the most advanced DSP electronics and auto-alignment mechanism. The FVS-6000 not only allows you to easily obtain fingerprint VCD spectra, but also has several unique features such as a measurement range extension option of $4000 \sim 750 \text{ cm}^{-1}$. Since the CD signals in the infrared region are one or more orders of magnitude lower than ECD signals in the UV-Vis region, high sensitivity and stability are required for a VCD spectrometer. The FVS-6000 is the VCD spectrometer of choice for highly sensitive VCD measurements.



Sensitivity and Stability

● Optimized optics

The FVS-6000 employs a 28° interferometer to reduce polarization effects and reflective optics to eliminate artifacts caused by birefringence for high-quality VCD data. The isolated mount of the IR light source and the thermal control of the PEM crystal offers long-term baseline stability.

● VCD auto alignment

The auto alignment function optimizes the VCD optics to reduce linear anisotropy artifacts. No manual alignment is required when the detector and/or optical filters are exchanged.

● Lock-in detection by DSP

Phase-sensitive lock-in detection by DSP (Digital Signal Processing) with a sophisticated algorithm optimized for VCD offers drastic improvement of the S/N ratio.

● Purge capability

The optics, sample chamber and detector housing are all purgeable to obtain high-precision data free of environmental interference.

Features

● Wide band coverage

The FVS-6000 can provide spectra over a wide spectral range from $3200 \sim 850 \text{ cm}^{-1}$ as standard. Optional detectors and filters can extend the range to $4000 \sim 750 \text{ cm}^{-1}$, enabling the measurement of OH and NH bands. Detectors can be easily exchanged with no manual alignment, and optical filters can be switched by a PC-controlled 6-position filter wheel.

● Narrow band mode

The narrow band mode allows measurements of small peaks at targeted absorption bands with high sensitivity by using optional band filters.

● Variety of sample cells


A thermostatted cell holder, 3-position sample shuttle, and range of small-volume liquid cells are available as options.

● Spectra Manager™ II

Spectra Manager II software, the cross-platform control and analysis package for all JASCO spectrometers, offers quick and easy VCD data acquisition and analysis. Theoretically calculated VCD intensities from a commercially available calculation software can be plotted as simulated spectra to be compared with the experimental data for further tertiary structure analysis of molecules.


Long-term stability

VCD and IR Spectra of pinene



Narrow band mode

VCD and IR Spectra of camphor



Specifications

Instrumentation	VCD instrument
Measurement range	Standard: 3200 ~ 850 cm ⁻¹ Option: 4000 ~ 750 cm ⁻¹
Resolution	16 ~ 0.5 cm ⁻¹
Noise level (measurement conditions)	less than 8 x 10 ⁻⁶ ΔAbs (4 cm ⁻¹ , 20 min accumulation)
Light source	High intensity ceramic source, insulated mount
Detector	Standard: MCT-V (3200 ~ 850 cm ⁻¹) Option: MCT-C (1000 ~ 750 cm ⁻¹) InSb (4000 ~ 2000 cm ⁻¹)
Lock-in detection	DSP
PC communication	USB2.0 for main unit RS-232C for PEM driver
Interferometer	28° Michelson interferometer
Optical filter	Standard: 3200 ~ 2000 cm ⁻¹ , 2000 ~ 850 cm ⁻¹ (6-position automatic switching) Option: 4000 ~ 2700 cm ⁻¹ , 1000 ~ 750 cm ⁻¹ , 1850 ~ 1550 cm ⁻¹
Polarizer	KRS-5
Options	3-position automated sample shuttle Fixed liquid cell, BaF ₂ or CaF ₂ window, 25 ~ 1000 μm cell path Demountable liquid cell, BaF ₂ or CaF ₂ window, including 50 and 100 μm cell spacer Temperature control cell (Peltier cooled)
Dimensions	590(W) x 670(D) x 313(H) mm
Weight	65 kg



*The VFT-4000 VCD attachment for the FT/IR-4000, 6000 Series FT-IR spectrometers is also available.

VFT-4000 VCD attachment with the FT/IR-6300 FT-IR spectrometer



● Specifications are subject to change without notice.

JASCO INTERNATIONAL CO., LTD.

4-21, Sennin-cho 2-chome, Hachioji, Tokyo 193-0835, Japan
Tel: +81-42-666-1322 Fax: +81-42-665-6512 <http://www.jascoint.co.jp/english/>
Australia, China, Hong Kong, India, Indonesia, Iran, Korea, Malaysia, New Zealand,
Pakistan, Philippines, Russia, Singapore, South Africa, Taiwan, Thailand

JASCO INCORPORATED

28600 Mary's Court, Easton, MD 21601, U.S.A
Tel: +1-800-333-5272 +1-410-822-1220 Fax: +1-410-822-7526 <http://www.jascoinc.com>
U.S.A., Canada, Costa Rica, Mexico, Puerto Rico, Argentina, Brazil, Chile, Colombia, Paraguay,
Peru, Uruguay, Guatemala

JASCO EUROPE s.r.l.

Via Luigi Cadorna 1, 23894 Cremella (Lc), Italy
Tel: +39-039-9215811 Fax: +39-039-9215835 <http://www.jasco-europe.com>
JASCO Deutschland www.jasco.de, JASCO UK www.jasco.co.uk,
JASCO France www.jascofrance.fr, JASCO Benelux www.jasco.nl,
JASCO Spain www.jasco-spain.com, JASCO Scandinavia www.jascoscandinavia.se
Italy, Germany, U.K., France, Netherlands, Belgium, Luxembourg, Spain, Sweden, Norway, Denmark,
Austria, Finland, Greece, Hungary, Poland, Portugal, Romania, Switzerland, Algeria, Cyprus,
Egypt, Israel, Jordan, Kuwait, Lebanon, Morocco, Saudi Arabia, Syria, Tunisia, Turkey,
U.A.E., Yemen

For more information, please contact: